

Chlamydomonas adapted GFP (CrGFP)**Plasmid pCrGFP****Sequence of the CrGFP**

In the following sequence, all amino acids which have been altered from the wildtype GFP from *Aequorea victoria* are printed in bold face. The nucleotides which have been exchanged for codon adaptation are also printed in bold face. The first amino acids in italic (up to ...EDPP) are encoded by the multiple cloning site (MCS) of the expression wildtype and are thus not identical with the wildtype amino acid sequence. The expression in *E. coli* or *Chlamydomonas* has not been affected by this, however.

There are only six nucleotides missing compared to the wildtype gene (at positions 52 to 57), with the sequence "ATG GCC", which code for methionine and alanine. These six nucleotides have been replaced by the nucleotides 1 to 57 in the sequence below due to ligation into the expression vector.

```

1 ATG ACC ATG ATT ACG CCA AGC TTG CAT GCC TGC AGG TCG ACT CTA GAG 48
    M   T   M   I   T   P   S   L   H   A   C   R   S   T   L   E

49 GAT CCC CCC AAG GGC GAG GAG CTG TTC ACC GGT GTG GTC CCC ATC CTG 96
    D   P   P   K   G   E   E   L   F   T   G   V   V   P   I   L

97 GTG GAG CTG GAC GGC GAC GTG AAC GGC CAC AAG TTC TCC GTC TCC GGC 144
    V   E   L   D   G   D   V   N   G   H   K   F   S   V   S   G

145 GAG GGT GAG GGT GAC GCC ACC TAC GGC AAG CTG ACC CTG AAG TTC ATC 192
    E   G   E   G   D   A   T   Y   G   K   L   T   L   K   F   I

193 TGC ACC ACC GGC AAG CTG CCC GTG CCC TGG CCC ACC CTG GTC ACC ACC 240
    C   T   T   G   K   L   P   V   P   W   P   T   L   V   T   T

241 CTG ACC TAC GGT GTG CAG TGC TTC TCC CGC TAC CCC GAC CAC ATG AAG 288
    L   T   Y   G   V   Q   C   F   S   R   Y   P   D   H   M   K

289 CAG CAC GAC TTC TTC AAG TCC GCC ATG CCC GAG GGC TAC GTG CAG GAG 336
    Q   H   D   F   F   K   S   A   M   P   E   G   Y   V   Q   E

337 CGC ACC ATC TTC TTC AAG GAC GAC GGC AAC TAC AAG ACC CGC GCC GAG 384
    R   T   I   F   F   K   D   D   G   N   Y   K   T   R   A   E

385 GTC AAG TTC GAG GGC GAC ACC CTG GTG AAC CGC ATC GAG CTG AAG GGC 432
    V   K   F   E   G   D   T   L   V   N   R   I   E   L   K   G

433 ATC GAC TTC AAG GAG GAC GGC AAC ATC CTG GGC CAC AAG CTG GAG TAC 480
    I   D   F   K   E   D   G   N   I   L   G   H   K   L   E   Y

481 AAC TAC AAC TCC CAC AAC GTG TAC ATC ATG GCC GAC AAG CAG AAG AAC 528
    N   Y   N   S   H   N   V   Y   I   M   A   D   K   Q   K   N

529 GGC ATC AAG GTG AAC TTC AAG ATC CGC CAC AAC ATC GAG GAC GGC TCC 576
    G   I   K   V   N   F   K   I   R   H   N   I   E   D   G   S

577 GTG CAG CTG GCC GAC CAC TAC CAG CAG AAC ACC CCC ATC GGC GAT GGC 624
    V   Q   L   A   D   H   Y   Q   Q   N   T   P   I   G   D   G

625 CCC GTG CTG CTG CCC GAC AAC CAC TAC CTG TCC ATC CAG TCC GCC CTG 672
    P   V   L   L   P   D   N   H   Y   L   S   I   Q   S   A   L

```



673 TCC **AAG** GAC CCC **AAC** **GAG** AAG **CGC** GAC CAC ATG GTC **CTG** **CTG** GAG TTC 720
 S K D P N E K R D H M V L L E F
 721 GTC ACC GCT GCC GGC ATC ACC CAC GGC ATG GAC GAG CTG TAC AAG TAA 768
 V T A A G I T H G M D E L Y K *
 After the stop codon follows a 3' untranslated region (3'-UTR) with the sequence:
 769 GGA TCC CCG CTC CGT GTA AAT GGA GGC GCT CGT TGA TCT GAG CCT TGC 816
 817 CCC CTG ACG AAC GGC GGT GGA TGG AAG ATA CTG CTC TCA AAG TGC TGA 864
 865 AGC GGT AGC TTA GCT CCC CGT TTC GTG CTG ATC AGT CTT TTT CAA CAC 912
 913 GTA AAA AGC GGA GGA GTT TTG CAA TTT TGT TGG TTG TAA CGA TCC TCC 960
 961 GTT GAT TTT GGC CTC TTT CTC CAT GGG CGG GCT GGG CGT ATT GAA GCG 1008
 1009 GGT ACC GAG CTC GAA TTC 1026

Base composition of the coding region:


A	23,05 %	GC	61,33 %
C	34,64 %	AT	38,67 %
G	26,69 %		
T	15,63 %		

Amino acid exchanges from the wildtype gene from *Aequorea victoria*

Two amino acids have been exchanged in the gene, F64L and S65T. The exchanges have effects on three properties of the protein: The fluorescence intensity at an excitation wavelength of 480 nm is increased 30fold; the formation of the chromophore is accelerated 10 to 20fold; and the protein shows only one excitation maximum at 489 nm versus two maxima (at 395 nm and 470 nm) in the wildtype. The reduction to one maximum suites the protein for confocal microscopy, because most laser-scanning confocal microscopes (LSCM) work with lasers of 488 nm emission wavelength.

The plasmid pCrGFP

The algal CrGFP is carried by the plasmid pCrGFP. This plasmid is based on the vector pUC19. The CrGFP and a 3' untranslated region have been ligated into the MCS of pUC19. This region is the 3'-UTR of the rbcS2 gene from *Chlamydomonas reinhardtii*, which codes for a sub unit of the RUBISCO.





The 5' part of the MCS lies immediately downstream of P_{lac} . As the detailed map of the MCS shows, genes can be ligated immediately in front of the CrGFP by using the existing restriction sites, in order to verify their expression using the GFP fluorescence.

Map of the restriction sites of the vector

Enzymes that do not cut in pCrGFP:

AfIII	Apal	Ascl	AvrII	Bael	Bael	BbsI
BgIII	Bpu10I	BsaBI	BsmI	BspEI	BssHII	BstZ17I
Bsu36I	ClaI	DrdII	EagI	Eco47III	EcoRV	Fsel
HpaI	MluI	MscI	MunI	NheI	NotI	NruI
NsiI	NspV	PacI	PmeI	PmlI	PshAI	RsrII
SacII	SanDI	SexAI	SfiI	SgfI	SmaI	SmI
SnaBII	SpeI	SrfI	Sse8647I	StuI	SunI	Tth111I
UbaDI	UbaEI	Xhol	XmaI			

Enzymes that cut, with position of restriction site:

AatII	1	Agel	1	AhdI	1	AlwNI	1
BcgI	1	BcgI	1	BclI	1	BplI	1
Bpu1102I	1	Bsal	1	BspGI	1	BspLU11I	1
BspMI	1	BstAPI	1	BstEII	1	BstXI	1
DraIII	1	EcoNI	1	EcoRI	1	HindIII	1
KpnI	1	NarI	1	Ncol	1	NdeI	1
NgoAIV	1	Pfl1108I	1	PstI	1	RleAI	1
SacI	1	Sall	1	SapI	1	SbfI	1
Scal	1	SgrAI	1	SphI	1	SspI	1
XbaI	1	XcmI	1	XmnI	1	AclI	2
BamHI	2	BciVI	2	Bfil	2	BglI	2
Bsbl	2	BseRI	2	BsrDI	2	BsrGI	2
FspI	2	HgiEII	2	PflMI	2	PvuI	2
ApaLI	3	BsaXI	3	BsmBI	3	DraI	3
DrdI	3	EarI	3	Pvull	3	RcaI	3
TaqII	3	TaqII	3	VspI	3	Bce83I	4
BpmI	4	BsgI	4	Bsp24I	4	Bsp24I	4
BsrBI	4	BssSI	4	Eco57I	5	EciI	6
AcellI	7						


Enzyme	Position
AatII	70
Agel	2294
AhdI	989
AlwNI	1468
BcgI	449
BcgI	485
BclI	3111
BplI	2787
Bpu1102I	3091
Bsal	923
BspGI	2323
BspLU11I	1877
BspMI	2245
BstAPI	3457
BstEII	2449
BstXI	2314
DraIII	2716
EcoNI	2891

Enzyme	Position						
EcoRI	3239						
HindIII	2236						
KpnI	3231						
NarI	3401						
Ncol	3199						
NdeI	3453						
NgoAIV	2948						
Pfl1108I	970						
PstI	2252						
RleAI	2726						
SacI	3237						
Sall	2254						
SapI	1997						
SbfI	2252						
Scal	508						
SgrAI	2294						
SphI	2246						
SspI	184						
XbaI	2260						
XcmI	3206						
XmnI	389						
AclI	387	760					
BamHI	2266	2987					
BciVI	152	1679					
BfiI	948	3280					
BglII	871	3391					
Bsbl	3127	3519					
BseRI	2297	3156					
BsrDI	755	927					
BsrGI	2718	2976					
FspI	766	3381					
HgiEII	1295	3453					
PflMI	2689	2965					
PvuI	620	3362					
ApaLI	317	1563	3458				
BsaXI	2023	2347	2704				
BsmBI	2359	3590	3636				
DraI	411	1103	1122				
DrdI	1775	2921	3544				
EarI	193	1997	3343				
PvuII	2057	2800	3331				
RcaI	44	149	1157				
TaqI	455	642	1981				
TaqII	287	302	2972				
VspI	814	2049	2108				
Bce83I	385	1253	1492	1792			
BpmI	920	2336	2711	2951			
BsgI	2393	2492	2567	2816			
Bsp24I	89	1178	1356	2850			
Bsp24I	52	1215	1393	2887			
BsrBI	147	1948	2189	2995			
BssSI	17	320	1704	3655			
Eco57I	323	1333	2420	2663	3100		
EcI	833	1661	1807	2527	2771	2884	
AclII	689	1929	2277	2310	2631	2964	3585



The multiple cloning site of pCrGFP

The CrGFP has been inserted in frame with the MCS, so that the protein can be expressed in *E. coli* immediately just by transformation.



The two *BamHI* sites allow the excision of the CrGFP without the 3'-UTR.

Transformation of *E. coli* with pCrGFP

To transform *E. coli* with pCrGFP in its original or a modified version, you need competent *E. coli* cells. You find detailed instructions to prepare such cells for electrical or chemical transformation in "Current Protocols in Molecular Biology" (Ausubel et al., 1990) or in "Molecular Biology: A Laboratory Manual" (Sambrook et al. 1989).

Map of the restriction sites of the CrGFP gene

This map contains all restriction sites which cut inside the CrGFP gene including the 3'-UTR. The gene starts at nucleotide 2218 of the plasmid and ends at nucleotide 2987. Then follows the 3'-UTR, which ends at nucleotide 3219.

Enzymes that do not cut in CrGFP:

AatII	AclI	AfIII	AhdI	AlwNI	Apal	ApaLI
Ascl	AvrII	Bael	Bael	BbsI	Bce83I	Bcgl
Bcgl	BciVI	Bfil	BgII	BgIII	Bpu10I	Bsal
BsAbI	BsmI	BspEI	BspLU11I	BsrDI	BssHII	BssSI
BstAPI	BstZ17I	Bsu36I	Clal	Dral	DrdII	EagI
EarI	Eco47III	EcoRV	Fsel	Fspl	HgiEII	Hpal
Mlul	Mscl	MunI	NarI	Ndel	Nhel	NotI
NruI	Nsil	NspV	PacI	Pfl1108I	PmeI	PmlI
PshAI	PvuI	RsrII	SacII	SanDI	SapI	Scal
SexAI	SfiI	SgfI	SmaI	Smil	SnaBI	Spel
SrfI	Sse8647I	SspI	StuI	SunI	TaqI	Tth111I
UbaDI	UbaEI	VspI	Xhol	XmaI	XmnI	

Enzymes that cut, with position of restriction sites:

The positions refer to the position of the gene in the plasmid. The restriction map below starts with nucleotide 1 of the gene (the A of the start codon), which lies at position 2218 in the plasmid.

AgeI	1	BclI	1	BplI	1	Bpu1102I	1
Bsbl	1	BsmBI	1	Bsp24I	1	Bsp24I	1
BspGI	1	BspMI	1	BsrBI	1	BstEII	1
BstXI	1	DraIII	1	DrdI	1	EcoNI	1
HindIII	1	Ncol	1	NgoAIV	1	PstI	1
Pvull	1	RleAI	1	Sall	1	SbfI	1
SgrAI	1	SphI	1	TaqII	1	XbaI	1
XcmI	1	BamHI	2	BsaXI	2	BseRI	2



Chlamydomonas adapted GFP

BsrGI	2	PflMI	2	BpmI	3	Ecil	3
Eco57I	3		AcellI	4	BsgI	4	

Enzyme	Position				
Agel	2294				
BcII	3111				
BpII	2787				
Bpu1102I	3091				
Bsbl	3127				
BsmBI	2359				
Bsp24I	2887				
Bsp24I	2850				
BspGI	2323				
BspMI	2245				
BsrBI	2995				
BstEII	2449				
BstXI	2314				
DraIII	2716				
DrdI	2921				
EcoNI	2891				
HindIII	2236				
Ncol	3199				
NgoAIV	2948				
PstI	2252				
Pvull	2800				
RleAI	2726				
Sall	2254				
SbfI	2252				
SgrAI	2294				
SphI	2246				
TaqII	2972				
XbaI	2260				
XcmI	3206				
BamHI	2266	2987			
BsaXI	2347	2704			
BseRI	2297	3156			
BsrGI	2718	2976			
PflMI	2689	2965			
BpmI	2336	2711	2951		
Ecil	2527	2771	2884		
Eco57I	2420	2663	3100		
AcellI	2277	2310	2631	2964	
BsgI	2393	2492	2567	2816	



Map of the CrGFP gene:


<p style="text-align: center;">Sall</p> <p>HindIII BspMI PstI \ \ \ 1 atgaccatgattacgcacgttgcacgcgtccaggtcgactcttagaggatcccccaag 60 tactggtaactaatgcgggtcgaacgtacggacgtccagctgagatctcttagggggttc ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>AgeI SgrAI AceIII BseRI BstXI BspGI BpmI \ \ \ \ \ \ 61 ggcgaggagctgttaccgggtgtggtccccatctgtctggagctggacggcgacgtgaac 120 ccggctctcgacaagggtggcacaccagggttaggacgacctcgacactgcccgtgcacttg ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>BsaXI BsmBI BsgI \ \ \ 121 ggccacaagttctccgtctccggcgagggtggatgggtacgcgccacctacggcaagctgacc 180 ccgggttcaagaggcagaggccgtccactccactcggtggatgccgttcgactgg ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>Eco57I BstEII \ \ 181 ctgaagttcatctgcaccacccggcaagctgcccgtgcccctggccaccctggtaccacc 240 gacttcaagtagacgtggccgttcacgggcacggggaccgggtgggaccagtggtag ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>BsgI \ 241 ctgacctacgggtgcagtgtttctccgtacccggaccatagaaggcggcacttc gactggatgcacacgtcacgaagaggcgatgggtgtacttcgtgtcaag ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>EciI BsgI \ \ 301 ttcaagtccggcatgcccggagggtacgtgcaggaggcgaccatcttcaaggacgac 360 aagttcaggccgtacgggtcccgatgcacgtcctcggtttagaaagaagttctgtctg ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>AceIII \ 361 ggcaactacaagacccggccggagggtcaagttcgaggggcgacaccctggtaaccgcac ccgtttagttctggcgccgtcccgatgtttccgtgtggaccacttggcgtag ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>Eco57I PflMI \ \ 421 gagctgaaggcatcgacttcaaggaggacggcaacatctggccacaagctggagtac ctgcacttccgttagctgaagttccctctgtccgtttagggaccgttgcacccat ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>BsrGI BsaXI BpmI DraIII RleAI \ \ \ \ \ 481 aactacaactccacaacgtgtacatcatggccgacaaggcagaagaacggcatcaagg ttgtatgttgggttgcacatgttagtaccggctgttgcgttccgttagtccac ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>EciI BplI PvuII BsgI \ \ \ \ 541 aacttcaagatccggccacaacatcgaggacggctccgtcgactggccgaccactacc ttgaagttctaggccgttgcacatgttagtaccggctgttgcgttccgttagtccac ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p> <p>Bsp24I \ 601 cagaacaccccatcgccgtggccctgtgtgtccgcacaaccactaccgtccatc gtcttgtgggttagccgtaccggggcagacgacgggtgttggatggacaggtag ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ * ^ *</p>	60 120 180 240 300 360 420 480 540 600 660
---	--






Fluorescence spectra

Emission spectrum with excitation at 488 nm



Excitation spectrum with emission at 512 nm






Microphotographs


Coloured versions of the photographs shown here can be found on our website under www.entelechon.com.



This photo shows cells of the single-celled alga *Clamydomonas reinhardtii* as seen under the light microscope. These cells have been transformed with a construct of the algal GFP and a resistance gene. The resistance gene binds to DNA, so that one can expect to find the GFP-resistance gene hybrid to be enriched in the nucleus.



Here, you can see the fluorescence of the cells shown above in the transmission at an excitation of 395 nm, with a bandpass emission filter of 510-520 nm. The emission filter is used to block the red chlorophyll fluorescence. You can clearly see the green fluorescence of the GFP, localized in a spherical, central region inside the cells.



In this picture, you see the fluorescence at an excitation of 330 nm, using a 520 nm longpass emission filter. The nucleus of the cells has been stained with ethidium bromide. There is red fluorescence of the chlorophyll from the chloroplast, as well as the yellowish fluorescence of ethidium bromide. The chloroplast has the typical cup-shaped outline, whereas the nucleus lies inside the chloroplast cup. The localization of GFP and ethidium bromide matches perfectly, showing that the GFP-resistance gene construct is located in the nucleus.

Technical support

If you have questions concerning our products or need help, please consult our technical support staff. We are also happy about any suggestions which help us to offer you better products and service. Contact us...

...via email:

support@entelechon.de

...via telephone:

++49-941-94 683 - 60

...via fax:

++49-941-94 683 - 66

...via mail:

Entelechon GmbH
Josef-Engert-Straße 9
93053 Regensburg
Germany

In general, the most convenient way to ask your questions is via email. This allows our support staff to examine your problem in detail and then call you back. You don't lose time with complicated explanations on the telephone, but instead get fast solutions to your problems.

License

By using the plasmid **pCrGFP** (product) delivered by Entelechon GmbH (licenser) or any products developed from pCrGFP, the user (licensee) agrees with this license.

Entelechon GmbH grants to the licensee the right without temporal limitation to use the plasmid **pCrGFP**.

The licensee is forbidden to do any of the following:

- Passing the product or its derivatives on to a third party, with or without payment.
- Production of illegal products or organisms, or products or organisms for an illegal application with the help of this product or its derivatives, or with the help of information from this product or its derivatives, or with the help of information supplied by the licensor with the product.
- Passing the information supplied with the product on to a third party, with or without payment.

In terms of this license, derivatives of the product are:



- Biomolecules, which have been gained from reproducing the product
- Biomolecules, which are a part of the product (e.g. partial sequences)
- Biomolecules, which have been gained from combining the product with other biomolecules.
- Biomolecules which are derived from this product by transcription and/or translation or other processes, which allow to create a related sequence from the specified sequence of a biomolecule.

Limited Warranty

Entelechon GmbH is committed to providing our customers with high-quality products and services. Our goal is to ensure that our customers are 100% satisfied with our products and services.

Entelechon GmbH warrants that all of its products will perform according to the specification in the product information supplied with the products. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Entelechon GmbH's liability only to the cost of the product. No warranty is granted for products beyond their expiration date. No warranty is applicable unless all product components are stored and handled in accordance with instructions. Entelechon reserves the right to select the method(s) used to analyze a product in case of reported errors.

Entelechon GmbH makes every effort to ensure the accuracy of its publications, but realizes that occasional typographical or other errors in published texts, descriptions, files and computer programs can not be excluded. Therefore, Entelechon GmbH makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our support staff (support@entelechon.de).

Entelechon GmbH assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Entelechon GmbH
Josef-Engert-Straße 9
93053 Regensburg
Tel. (0941) 94 683 - 60
Fax. (0941) 94 683 - 66
web: www.entelechon.com
email: support@entelechon.de